

RHBDF1 Rabbit pAb

货号: **AYP21085**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	97kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Mouse heart,Mouse kidney,Rat lung,Rat heart,Rat kidney
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 200-400 of human RHBDF1 (NP_071895.3).
序列	HRLPRRRKRESVAKMSFRAAAAALMKGRSVRDGTFRRAQRSSFTPASFLEEDTTDFPDELDTSSFFAREGILHEELSTYPDEV FESPSEAAALKDWEKAPEQADLTGGALDRSELSHLMPLERGWKQKEGAAAPQPKVRLRQEVVSTAGPRRGQRIAVP VRKLFAREKRPYGLGMVGRILTNRITYRKRIDSFVKRQIEDMD

靶点信息

研究背景	Regulates ADAM17 protease, a sheddase of the epidermal growth factor (EGF receptor ligands and TNF, t hereby plays a role in sleep, cell survival, proliferation, migration and inflammation. Does not exhibit any protease activity on its own.
基因ID	64285

基因名	RHBDF1
Swiss	Q96CC6
别名	C16orf8;Dist1;EGFR-RS;gene-89;gene-90;hDist1;RHBDF1

产品验证



Western blot analysis of RHBDF1 expressed in Mouse lung,Rat lung,Mouse heart using RHBDF1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 12 0s.

实验步骤

访问官网浏览详情: www.ablybio.cn