

NSRP1 Rabbit pAb

货号: **AYP21090**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	75kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,Jurkat,Mouse lung,Rat brain
细胞定位	nuclear speck,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human NSRP1 (NP_115517.1).
序列	MAIPGRQYGLILPKKTQQLHPVLQKPSVFGNDSDDDDDETSVSESLQREAAKKQAMKQTKLEIQKALAEDATVYEYDSIYDE MQKKKEENPKLLLGGKDRKPKYIHNLLKAVEIRKKEQEKRMEKKIQREREMEKGEFDDKEAFVTSAYKKKLQERAEEREK RAAALEACLDVTKQKDLSGFYRHLLNQAVGEEVPKCSFREARSGIKEEKSRGFSNEVSSKNRIPQEKCILQTDVKVEENP DADSDFDKSSADDEIEETRVNCRREKVIETPENDF

靶点信息

研究背景	RNA-binding protein that mediates pre-mRNA alternative splicing regulation.
基因ID	84081
基因名	NSRP1

Swiss	Q9H0G5
别名	CCDC55;HSPC095;NSRp70;NSRP1

产品验证



Western blot analysis of NSRP1 expressed in Mouse lung,Rat brain,HeLa using NSRP1 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn