

— ABLYBIO, Help Your Research



NSRP1 Rabbit pAb

货号: **AYP21090**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:2000 |
| 理论分子量 | |
| 实测分子量 | 75kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | HeLa,293T,Jurkat,Mouse lung,Rat brain |
| 细胞定位 | nuclear speck,nucleoplasm,nucleus |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|---|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human NSRP1 (NP_115517.1). |
|------|---|

靶点信息

| | |
|------|---|
| 研究背景 | RNA-binding protein that mediates pre-mRNA alternative splicing regulation. |
|------|---|

| | |
|-------|--|
| 基因ID | 84081 |
| 基因名 | NSRP1 |
| Swiss | Q9H0G5 (https://www.uniprot.org/uniprotkb/Q9H0G5/entry) |
| 别名 | CCDC55,HSPC095,NSrp70,NSRP1,NSRP1 Rabbit pAb,Coiled-coil domain-containing protein 55,Nuclear speckle-related protein 70 |

产品验证



Western blot analysis of NSRP1 expressed in Mouse lung,Rat brain,HeLa using NSRP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)