

# TOX2 Rabbit pAb

货号: **AYP21096**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	<a href="#">WB</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,U-251MG,Rat lung,Rat thymus
细胞定位	nucleoplasm,nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human TOX2 ( NP_116272.1).
序列	MSDGNPELLSTSQTYNGQSENNEDYEIPPITPPNLPEPSLLHLGDHEASYHSLCHGLTPNGLLPAYSYQAMDLP AIMVSNM LAQDSHLLSGQLPTIQEMVHSEVAAYDSGRPGLLGRPAMLASHMSALSQSQLISQMGI

## 靶点信息

研究背景	Putative transcriptional activator involved in the hypothalamo-pituitary-gonadal system.
基因ID	84969
基因名	TOX2

<b>Swiss</b>	Q96NM4
别名	C20orf100;GCX-1;GCX1;dj1108D11.2;dj495O3.1;TOX2

## 产品验证



Western blot analysis of TOX2 expressed in Rat lung,U-87MG,U-251MG using TOX2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)