

— ABLYBIO, Help Your Research



MED19 Rabbit pAb

货号: **AYP21140**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	20kDa/26kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Mouse liver
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 60-160 of human MED19 (NP_703151.1).
------	--

靶点信息

研究背景	The protein encoded by this gene is a subunit of the Mediator complex, which binds to gene-specific regulatory factors and provides support for the basal RNA polymerase II transcription machinery. This gene has been implicated in the growth of several types of cancer, and inhibition of its expression inhibits the growth and spread of these cancers. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2015]
基因ID	219541
基因名	MED19
Swiss	A0JLT2 (https://www.uniprot.org/uniprotkb/A0JLT2/entry)
别名	DT2P1G7,LCMR1,MED19AS,MED19,MED19 Rabbit pAb,Lung cancer metastasis-related protein 1,Mediator complex subunit 19

产品验证



Western blot analysis of MED19 expressed in U-87MG, Mouse liver using MED19 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)