

STAT5A/STAT5B Rabbit pAb

货号: AYP21159

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	90kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	mouse kidney
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 700-787 of human STAT5A/STAT5B (NP_036580.2).
序列	VKPQIKQVVPEFVNASADAGGGSATYMDQAPSPAVCPQAHYVMYPQNPDSVLDTGDFDLEDTMDVARRVEELLGRP MDSQWIPHAQS

靶点信息

研究背景	Carries out a dual signal transduction and activation of transcription. Mediates cellular responses to the cytokine KITLG/SCF and other growth factors. Mediates cellular responses to ERBB4. May mediate cellular responses to activated FGFR1, FGFR2, FGFR3 and FGFR4. Binds to the GAS element and activates PRL-induced transcription. Regulates the expression of milk proteins during lactation.
基因ID	6776,6777

基因名	STAT5A,STAT5B
Swiss	
别名	MGF;STAT5;stat5a/stat5b

产品验证



Western blot analysis of STAT5A/STAT5B expressed in mouse kidney using STAT5A/STAT5B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn