

— ABLYBIO, Help Your Research



PEBP4 Rabbit pAb

货号: **AYP21171**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Human plasma,Rat testis
细胞定位	extracellular region
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 90-242 of mouse PEBP4 (NP_082836.2).
------	--

靶点信息

研究背景	Promotes AKT phosphorylation, suggesting a possible role in the PI3K-AKT signaling pathway.
基因ID	73523
基因名	Pebp4
Swiss	Q9D9G2 (https://www.uniprot.org/uniprotkb/Q9D9G2/entry)
别名	PEBP-4,PEBP4,PEBP4 Rabbit pAb

产品验证



Western blot analysis of PEBP4 expressed in Human plasma,Rat testis using PEBP4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)