

PEBP4 Rabbit pAb

货号: **AYP21171**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IHC IF/ICC |
| 推荐浓度 | WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 |
| 理论分子量 | |
| 实测分子量 | 26kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Human plasma,Rat testis |
| 细胞定位 | extracellular region |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|---|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 90-242 of mouse PEBP4 (NP_082836.2). |
| 序列 | RQKITAWQAPIVKFHTALDGALYLLVMVDPDAPSRSNPVMKYWRHWLVSNITGADMKSGSIRGNVLSDYSPPTPPPETGL HRYQFFVYLQGDRDISLSVEEKADLGGWNLDKFLQQYGLRDPDTSTQFMTQFDEELSSEFGRINDDQEQFNQK |

靶点信息

| | |
|------|---|
| 研究背景 | Promotes AKT phosphorylation, suggesting a possible role in the PI3K-AKT signaling pathway. |
| 基因ID | 73523 |
| 基因名 | |

| | |
|-------|--------------|
| Swiss | Q9D9G2 |
| 别名 | PEBP-4;PEBP4 |

产品验证



Western blot analysis of PEBP4 expressed in Human plasma,Rat testis using PEBP4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn