

MST1/STK4 Rabbit pAb

货号: AYP21273

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	52kDa/55kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 303-487 of human MST 1/STK4 (NP_006273.1).
序列	KRQESQQREVDQDDEENSEEDEMDSGTMVRAVGDGEMGTVRVASTMTDGANTMIEHDDTLPSQLGTMVINAEDDEEEEG TMKRRDET MQPAKPSFLEYFEQKEKENQINSFGKSVPGPLKNSSDWKIPQGDYEF LKSWTVEDLQKRLALDPMMEQE IEEIRQKYQSKRQPILDAIEAKRRRQQNF

靶点信息

研究背景	The protein encoded by this gene is a cytoplasmic kinase that is structurally similar to the yeast Ste20p kinase, which acts upstream of the stress-induced mitogen-activated protein kinase cascade. The encoded protein can phosphorylate myelin basic protein and undergoes autophosphorylation. A caspase-cleaved fragment of the encoded protein has been shown to be capable of phosphorylating histone H2B. The particular phosphorylation catalyzed by this protein has been correlated with apoptosis, and it's possible that this protein induces the chromatin condensation observed in this process.
基因ID	6789
基因名	STK4
Swiss	Q13043
别名	KRS2;MST1;YSK3;STK4

产品验证



Western blot analysis of MST1/STK4 expressed in HeLa using MST1/STK4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn