

— ABLYBIO, Help Your Research



ARHGAP4 Rabbit pAb

货号: **AYP21321**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	105kDa/109kDa
实测分子量	105kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat
细胞定位	cytoplasm,cytosol,Golgi apparatus,nucleoplasm
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 900 to the C-terminus of human ARHGAP4 (NP_001158213.1).
------	---

靶点信息

研究背景	This gene encodes a member of the rhoGAP family of proteins which play a role in the regulation of small GTP-binding proteins belonging to the RAS superfamily. The protein encoded by the orthologous gene in rat is localized to the Golgi complex and can redistribute to microtubules. The rat protein stimulates the activity of some Rho GTPases in vitro. Genomic deletions of this gene and a neighboring gene have been found in patients with nephrogenic diabetes insipidus. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	393
基因名	ARHGAP4
Swiss	P98171 (https://www.uniprot.org/uniprotkb/P98171/entry)
别名	C1,RGC1,RhoGAP4,SrGAP4,p115,ARHGAP4,ARHGAP4 Rabbit pAb,Rho-GAP hematopoietic protein C1,Rho-type GTPase-activating protein 4,KIAA0131

产品验证



Western blot analysis of ARHGAP4 expressed in Jurkat using ARHGAP4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)