

# OSR1 Rabbit pAb

货号: **AYP21352**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	29kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NCI-H460,Cal-27,A-549,22Rv1
细胞定位	cytosol,nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-266 of human OSR1 ( NP_660303.1).
序列	MGSKTLPAPVPIHPSLQLTNYSFLQAVNGLPTVPSDHLPLNYGFSALHAVHLHQWTLGYPAMHLPRSSFSKVPGTVSSLV DARFQLPAFPWFPHVIQPKPEITAGGSVPALKTKPRDFANLALAATQEDPAKLGRGEGPGSPAGGLGALLDVTKLSPEKK PTRGRLPSKTKKEFVCKFCGRHFTKSYNLLIHERHTHTDERPYTCDICHKAFRRQDHLRDHRYIHSKEKPFKQCQCGKGFQC SRTLAVHKTLSQVKELKTSKIKC

## 靶点信息

研究背景	Transcription factor that plays a role in the regulation of embryonic heart and urogenital development.
基因ID	130497

基因名	OSR1
Swiss	Q8TAX0
别名	ODD;OSR1

## 产品验证



Western blot analysis of OSR1 expressed in NCI-H460, Cal-27, A-549 using OSR1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)