

— ABLYBIO, Help Your Research



OSR1 Rabbit pAb

货号: **AYP21352**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IHC |
| 推荐浓度 | WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 |
| 理论分子量 | 29kDa |
| 实测分子量 | 30kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | NCI-H460,Cal-27,A-549,22Rv1 |
| 细胞定位 | cytosol,nucleus |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-266 of human OSR1 (NP_660303.1). |
|------|--|

靶点信息

| | |
|------|---|
| 研究背景 | Transcription factor that plays a role in the regulation of embryonic heart and urogenital development. |
|------|---|

| | |
|-------|--|
| 基因ID | 130497 |
| 基因名 | OSR1 |
| Swiss | Q8TAX0 (https://www.uniprot.org/uniprotkb/Q8TAX0/entry) |
| 别名 | ODD,OSR1,OSR1 Rabbit pAb |

产品验证



Western blot analysis of OSR1 expressed in NCI-H460, Cal-27, A-549 using OSR1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)