

OSR2 Rabbit pAb

货号: **AYP21363**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,293T,Mouse testis,Mouse pancreas,Rat large intestine,Rat kidney
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 250-312 of human OSR2 (NP_001135934.1).
------	--

靶点信息

研究背景	OSR2 is a mammalian homolog of the Drosophila odd-skipped family of transcription factors (Lan et al., 2004 [PubMed 15175245]).[supplied by OMIM, Mar 2008]
------	---

基因ID	116039
基因名	OSR2
Swiss	Q8N2R0 (https://www.uniprot.org/uniprotkb/Q8N2R0/entry)
别名	OSR2,OSR2 Rabbit pAb

产品验证



Western blot analysis of OSR2 expressed in Mouse testis,Rat large intestine,SKOV3 using OSR2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)