

ANXA3 Rabbit pAb

货号: **AYP21483**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	36kDa
实测分子量	34KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	axon,cytoplasm,dendrite,extracellular exosome,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human ANXA3 (NP_005130.1).
序列	MASIWVGHRTVRDYPDFSPVDAEAIQKAIRGIGTDEKMLISILTERSNAQRQLVKEYQAAYGKELKDDLKGDLSGHFE HLMVALVTPPAVFDKQLKSMKGAGTNEALIEILTRTSRQMKDISQAYYTVYKSLGDDISSETSGDFRKALLTLADG RRDESLKVDEHLAKQDAQILYKAGENRWGTDEDKFTIELCLRSFPQLKLTDFEYRNISQKDIVDSIKGELSGHFEDLLLAIVN CVRNTPAFLAERLHRALKGIGTDEFTLNRIMVSRSEIDLLDIRTEFKKHGYSLYSAIKSDTSGDYEITLLKICGGDD

靶点信息

研究背景	This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions in the inhibition of phospholipase A2 and cleavage of inositol 1,2-cyclic phosphate to form inositol 1-phosphate. This protein may also play a role in anti-coagulation.
------	---

基因ID	306
基因名	ANXA3
Swiss	P12429
别名	ANXA3;ANX3

产品验证



Western blot analysis of ANXA3 expressed in HeLa using ANXA3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn