

— ABLYBIO, Help Your Research



# TXLNA Rabbit pAb

货号: **AYP21499**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	61kDa
实测分子量	80KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	cytoplasm,cytosol,extracellular region
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TXLNA (NP_787048.1).
------	---

## 靶点信息

研究背景	May be involved in intracellular vesicle traffic and potentially in calcium-dependent exocytosis in neuroendocrine cells.
------	---

基因ID	200081
基因名	TXLNA
Swiss	P40222 ( <a href="https://www.uniprot.org/uniprotkb/P40222/entry">https://www.uniprot.org/uniprotkb/P40222/entry</a> )
别名	IL14, TXLN, TXLNA, TXLNA Rabbit pAb

## 产品验证



Western blot analysis of TXLNA expressed in HeLa using TXLNA Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)