

TXLNA Rabbit pAb

货号: **AYP21499**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	61kDa
实测分子量	80KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	cytoplasm,cytosol,extracellular region
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TXLNA (NP_787048.1).
序列	MKNQDKKNGAAKQSNPKSSPGQPEAGPEGAQERPSQAAPAVEAEGPGSSQAPRKPEGAQARTAQSGALRDVSEELSR QLEDILSTYCVDDNNQGGPGEDGAQGEPAEPEDAEKSRTYVARNGEPEPTPVVNGEKEPSKGDPNTEEIRQSDEVGDRDH RRPQ

靶点信息

研究背景	May be involved in intracellular vesicle traffic and potentially in calcium-dependent exocytosis in neuroendocrine cells.
基因ID	200081

基因名	TXLNA
Swiss	P40222
别名	IL14;TXLN;TXLNA

产品验证



Western blot analysis of TXLNA expressed in HeLa using TXLNA Rabbit pAb at 1:1000. Secondary antibody : HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn