

# TRMT2A/HTF9C Rabbit pAb

货号: **AYP21583**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:1000 <b>IHC:</b> 1:50 - 1:200
理论分子量	61kDa/68kDa
实测分子量	80KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	cytosol
纯化	Affinity purification

## 抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human TRMT2A/HTF9C (NP_073564.3).
序列	LASQHLVAILDPPRAGLHASKVILAIRRAKNLRRLLYVSCNPRAAMGNFVDLCRAPSNRVKGIPIFRPVKAVAVDLFPQTPHCEMLILFERVEHPNGTGVLGPHSPPAQPTGPPDNTLQETGTFPSS

## 靶点信息

研究背景	The protein encoded by this gene is of unknown function. However, it is orthologous to the mouse Trmt2a gene and contains an RNA methyltransferase domain. Expression of this gene varies during the cell cycle, with aberrant expression being a possible biomarker in certain breast cancers. Several transcript variants encoding two different isoforms have been found for this gene.
基因ID	27037

基因名	TRMT2A
Swiss	Q8IZ69
别名	TRMT2A;HTF9C

## 产品验证



Western blot analysis of TRMT2A/HTF9C expressed in HeLa using TRMT2A/HTF9C Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)