

— ABLYBIO, Help Your Research



ITSN1 Rabbit pAb

货号: **AYP21613**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	137kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG, SKOV3
细胞定位	clathrin-coated pit, cytoplasm, cytosol, neuron projection, nuclear envelope, plasma membrane, presynapse, recycling endosome
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 750-850 of human ITSN1 (NP_001317941.1).
------	---

靶点信息

研究背景	
------	--

基因ID	6453
基因名	ITSN1
Swiss	Q15811 (https://www.uniprot.org/uniprotkb/Q15811/entry)
别名	ITSN,SH3D1A,SH3P17,ITSN1,ITSN1 Rabbit pAb,SH3 domain-containing protein 1A

产品验证



Western blot analysis of ITSN1 expressed in U-251MG,SKOV3 using ITSN1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)