

— ABLYBIO, Help Your Research



# ZNF645 Rabbit pAb

货号: **AYP21665**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	
实测分子量	49kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Rat testis
细胞定位	cytoplasm
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 346-425 of human ZNF645 (NP_689790.1).
------	--

## 靶点信息

研究背景	
------	--

基因ID	158506
基因名	ZNF645
Swiss	Q8N7E2 ( <a href="https://www.uniprot.org/uniprotkb/Q8N7E2/entry">https://www.uniprot.org/uniprotkb/Q8N7E2/entry</a> )
别名	CT138,HAKAIL,ZNF645,ZNF645 Rabbit pAb,Cbl proto-oncogene-like protein 2,RING-type E3 ubiquitin transferase ZNF645,Zinc finger protein 645,c-Cbl-like protein 2

## 产品验证



Western blot analysis of ZNF645 expressed in Mouse testis,Rat testis using ZNF645 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)