

— ABLYBIO, Help Your Research



DTX3L Rabbit pAb

货号: **AYP21666**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	84kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse thymus
细胞定位	cytoplasm,cytosol,lysosome,nucleoplasm,nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 601-740 of human DTX3 L (NP_612144.1).
------	--

靶点信息

研究背景	DTX3L functions as an E3 ubiquitin ligase (Takeyama et al., 2003 [PubMed 12670957]).[supplied by OMI M, Nov 2009]
基因ID	151636
基因名	DTX3L
Swiss	Q8TDB6 (https://www.uniprot.org/uniprotkb/Q8TDB6/entry)
别名	BBAP,DTX3L,DTX3L Rabbit pAb,B-lymphoma- and BAL-associated protein,Protein deltex-3-like,RING-type E3 ubiquitin transferase DTX3L,Rhysin-2

产品验证



Western blot analysis of DTX3L expressed in Mouse thymus using DTX3L Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)