

# PPP2R1A/PPP2R2A/PPP2R1B Rabbit pAb

货号: **AYP21671**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	
实测分子量	55kDa/66kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2,K-562,NIH/3T3,C6,Mouse liver,Mouse lung,Rat brain,Rat testis
细胞定位	cytoplasm,cytosol,dendrite,extracellular exosome,glutamatergic synapse,lateral plasma membrane,micro tubule cytoskeleton,mitochondrion,neuron projection,nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-161 of human PPP2R1 A/PPP2R2A/PPP2R1BPPP2R1A (NP_055040.2).
序列	MAAADGDDSLYPIAVLIDELRNEDVQLRLNSIKKLSTIALALGVERTRSELLPFLTDTIYDEDEVLLALAEQLGTFTTLVGGPE YVHCLLPPELATVEETVVRDKAVESLRAISHEHSPSDLEAHFVPLVKRLAGGDWFTSRTSACGLFSVCYPRVSSA

## 靶点信息

研究背景	
基因ID	5518, 5520, 5519
基因名	PPP2R1A

<b>Swiss</b>	P30153,P63151,P30154
别名	PPP2R1A/PPP2R2A/PPP2R1B

## 产品验证



Western blot analysis of PPP2R1A/PPP2R2A/PPP2R1B expressed in Mouse liver,Rat brain,HeLa using PPP2R1A/PPP2R2A/PPP2R1B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TB ST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)