

— ABLYBIO, Help Your Research



PRELP Rabbit pAb

货号: **AYP21673**

产品信息

| | |
|-------|--|
| 反应 | Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:2000 |
| 理论分子量 | |
| 实测分子量 | 44KDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Mouse lung,Rat heart |
| 细胞定位 | extracellular exosome,extracellular region,extracellular space,extracellular vesicle,Golgi lumen,lysosomal lumen |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|---|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 200-360 of human PRELP (NP_002716.1). |
|------|---|

靶点信息

| | |
|------|--|
| 研究背景 | |
|------|--|

| | |
|-------|--|
| 基因ID | 5549 |
| 基因名 | PRELP |
| Swiss | P51888 (https://www.uniprot.org/uniprotkb/P51888/entry) |
| 别名 | MST161,MSTP161,SLRR2A,PRELP,PRELP Rabbit pAb,Proline-arginine-rich end leucine-rich repeat protein |

产品验证



Western blot analysis of PRELP expressed in Mouse lung,Rat heart using PRELP Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)