

— ABLYBIO, Help Your Research



SPAG11A Rabbit pAb

货号: **AYP22034**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	8kDa/9kDa/12kDa/13kDa
实测分子量	17kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat testis
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-50 of human SPAG11A (NP_47810 8.2).
------	--

靶点信息

研究背景	May have antimicrobial activity. May also play a role sperm maturation, storage, and protection.
------	--

基因ID	653423
基因名	SPAG11A
Swiss	Q6PDA7 (https://www.uniprot.org/uniprotkb/Q6PDA7/entry)
别名	SPAG11A,EDDM2A,HE2,sperm-associated antigen 11A,SPAG11A Rabbit pAb,Antimicrobial-like protein Bin-1b,Human epididymis-specific protein 2,Protein EP2,Sperm antigen HE2,BIN1B

产品验证



Western blot analysis of SPAG11A expressed in Rat testis using SPAG11A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)