

— ABLYBIO, Help Your Research



# Phospho-BRD4-S492/S494 Rabbit pAb

货号: AYP22350

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	DB WB
推荐浓度	<b>DB:</b> 1:500 - 1:1000 <b>WB:</b> 1:500 - 1:2000
理论分子量	80kDa/88kDa/152kDa
实测分子量	80kDa/88kDa/152kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2 cells, A549 cells, MCF-7 cells, Ramos cells, Jurkat cells, MDA-MB-231 cells, T-47D cells
细胞定位	Chromosome,Nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	A synthetic phosphorylated peptide around S492 & S494 of human BRD4 (NP_490597.1).
------	--

## 靶点信息

研究背景	The protein encoded by this gene is homologous to the murine protein MCAP, which associates with chromosomes during mitosis, and to the human RING3 protein, a serine/threonine kinase. Each of these proteins contains two bromodomains, a conserved sequence motif which may be involved in chromatin targeting. This gene has been implicated as the chromosome 19 target of translocation t(15;19)(q13;p13.1), which defines an upper respiratory tract carcinoma in young people. Two alternatively spliced transcript variants have been described.
基因ID	23476
基因名	BRD4
Swiss	O60885 ( <a href="https://www.uniprot.org/uniprotkb/O60885/entry">https://www.uniprot.org/uniprotkb/O60885/entry</a> )
别名	BRD4,CAP,HUNK1,HUNKI,MCAP,bromodomain containing 4,Phospho-BRD4-S492/S494 Rabbit pAb,Protein HUNK1,Phospho-BRD4-S492/S494

## 产品验证



Western blot analysis of Phospho-BRD4-S492/S494 expressed in HepG2 cells, A549 cells, MCF-7 cells using Phospho-BRD4-S492/S494 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))