

— ABLYBIO, Help Your Research



SLFN12 Rabbit pAb

货号: **AYP22372**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	67kDa
实测分子量	67KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y
细胞定位	cytosol, nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 320-390 of human SLFN 12 (NP_001275938.1).
------	--

靶点信息

研究背景	Predicted to act upstream of or within negative regulation of cell population proliferation.
------	--

基因ID	55106
基因名	SLFN12
Swiss	Q8IYM2 (https://www.uniprot.org/uniprotkb/Q8IYM2/entry)
别名	SLFN3,SLFN12 Rabbit pAb,SLFN12,Schlafen family member 12

产品验证



Western blot analysis of SLFN12 expressed in SH-SY5Y using SLFN12 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)