

— ABLYBIO, Help Your Research



RREB1 Rabbit pAb

货号: **AYP22378**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	181kDa
实测分子量	250KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,U-251MG,Rat lung
细胞定位	cytoplasm,extracellular exosome,fibrillar center,nuclear body,nuclear speck,nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 350-480 of human RREB 1 (NP_001003698.1).
------	---

靶点信息

研究背景	The protein encoded by this gene is a zinc finger transcription factor that binds to RAS-responsive elements (RREs) of gene promoters. It has been shown that the calcitonin gene promoter contains an RRE and that the encoded protein binds there and increases expression of calcitonin, which may be involved in Ras/Raf-mediated cell differentiation. Multiple transcript variants encoding several different isoforms have been found for this gene.
基因ID	6239
基因名	RREB1
Swiss	Q92766 (https://www.uniprot.org/uniprotkb/Q92766/entry)
别名	HNT,FINB,LZ321,Zep-1,RREB-1,RREB1 Rabbit pAb,RREB1,Finger protein in nuclear bodies,Raf-responsive zinc finger protein LZ321,Zinc finger motif enhancer-binding protein 1

产品验证



Western blot analysis of RREB1 expressed in Jurkat,U-251MG,Rat lung using RREB1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)