

— ABLYBIO, Help Your Research



Phospho-c-Jun-T91 Rabbit pAb

货号: **AYP22417**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	36KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,NIH/3T3
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	A synthetic phosphorylated peptide around T91 of human JUN (NP_002219.1).
------	---

靶点信息

研究背景	This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.
基因ID	3725
基因名	JUN
Swiss	P05412 (https://www.uniprot.org/uniprotkb/P05412/entry)
别名	AP-1,AP1,c-Jun,JUN,c,Phospho-c-Jun-T91 Rabbit pAb,Activator protein 1,Proto-oncogene c-Jun,Transcription factor AP-1 subunit Jun,V-jun avian sarcoma virus 17 oncogene homolog,p39

产品验证



Western blot analysis of Phospho-c-Jun-T91 expressed in 293T,NIH/3T3 using Phospho-c-Jun-T91 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)