

SARS-CoV-2 NSP10 Rabbit pAb

货号: AYP22632

产品信息

反应	SARS-CoV-2
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:2000 - 1:10000
理论分子量	
实测分子量	20-25KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-139 of coronavirus NS P10 (YP_009725306.1).
序列	AGNATEVPANSTVLSFCFAVDAAKAYKDYLASGGQPITNCVKMLCTHTGTGQAITVTPEANMDQESFGGASCCLYCRCH IDHPNPKGFCDLKGKVVQIPTTCANDPVGFTLKNTVCTVCGMWKGYGCSCDQLREPLQ

靶点信息

研究背景	[Replicase polyprotein 1ab]: Multifunctional protein involved in the transcription and replication of viral RNAs. Contains the proteinases responsible for the cleavages of the polyprotein.
基因ID	43740578
基因名	

Swiss	PODTD1
别名	

产品验证



Western blot analysis of SARS-CoV-2 NSP10 expressed in 293T using SARS-CoV-2 NSP10 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn