

— ABLYBIO, Help Your Research



CD96 Rabbit pAb

货号: **AYP22672**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	65kDa
实测分子量	100-200KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse spleen,Rat thymus
细胞定位	cytoplasm,plasma membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 270-370 of human CD96 (NP_937839.1).
------	--

靶点信息

研究背景	The protein encoded by this gene belongs to the immunoglobulin superfamily. It is a type I membrane protein. The protein may play a role in the adhesive interactions of activated T and NK cells during the late phase of the immune response. It may also function in antigen presentation. Alternative splicing generates multiple transcript variants encoding distinct isoforms.
基因ID	10225
基因名	CD96
Swiss	P40200 (https://www.uniprot.org/uniprotkb/P40200/entry)
别名	TACTILE,CD96 Rabbit pAb,CD96,Cell surface antigen CD96,T cell-activated increased late expression protein

产品验证



Western blot analysis of CD96 expressed in Jurkat, Mouse spleen, Rat thymus using CD96 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)