

IL13RA2 Rabbit pAb

货号: **AYP22718**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IHC IF/ICC |
| 推荐浓度 | WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 |
| 理论分子量 | 44kDa |
| 实测分子量 | 44kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | A-549,PC-3,Rat liver |
| 细胞定位 | Membrane,Single-pass type I membrane protein |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|---|
| 抗原信息 | A synthetic peptide corresponding to a sequence within amino acids 50-150 of human IL13RA2 (NP_000631.1). |
| 序列 | YLQWQPPLSLDHFKECTVEYELKYRNIGSETWKTIITKNLHYKDGFDLNKGIEAKIHTLLPWQCTNGSEVQSSWAETTYWIS PQGIPETKVQDMDCVYYNW |

靶点信息

| | |
|------|---|
| 研究背景 | The protein encoded by this gene is closely related to IL13RA1, a subunit of the interleukin 13 receptor complex. This protein binds IL13 with high affinity, but lacks cytoplasmic domain, and does not appear to function as a signal mediator. It is reported to play a role in the internalization of IL13. |
| 基因ID | 3598 |

| | |
|-------|------------------------------------|
| 基因名 | IL13RA2 |
| Swiss | Q14627 |
| 别名 | IL13RA2;CD213A2;CT19;IL-13R;IL13BP |

产品验证



Western blot analysis of IL13RA2 expressed in A-549,PC-3,Rat liver using IL13RA2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn