

— ABLYBIO, Help Your Research



HOXB6 Rabbit pAb

货号: **AYP22860**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	25kDa
实测分子量	25kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse lung
细胞定位	nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 30-110 of human HOXB 6 (NP_061825.2).
------	---

靶点信息

研究背景	This gene is a member of the Antp homeobox family and encodes a protein with a homeobox DNA-binding domain. It is included in a cluster of homeobox B genes located on chromosome 17. The encoded protein functions as a sequence-specific transcription factor that is involved in development, including that of lung and skin, and has been localized to both the nucleus and cytoplasm. Altered expression of this gene or a change in the subcellular localization of its protein is associated with some cases of acute myeloid leukemia and colorectal cancer.
基因ID	3216
基因名	HOXB6
Swiss	P17509 (https://www.uniprot.org/uniprotkb/P17509/entry)
别名	HOX2,HU-2,HOX2B,Hox-2.2,HOXB6 Rabbit pAb,HOXB6,Homeobox protein Hox-2.2,Homeobox protein Hox-2B,Homeobox protein Hu-2

产品验证



Western blot analysis of HOXB6 expressed in 293T, Mouse lung using HOXB6 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)