

— ABLYBIO, Help Your Research



SCRN1 Rabbit pAb

货号: **AYP22867**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IF/ICC |
| 推荐浓度 | WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200 |
| 理论分子量 | 46kDa |
| 实测分子量 | 50KDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | A-549,PC-3,U-251MG,Mouse brain,Rat brain |
| 细胞定位 | cytoplasm,nucleus,presynapse |
| 纯化 | Affinity purification |


抗原信息

| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 195-414 of human SCR N1 (NP_055581.3). |
|------|--|

靶点信息

| | |
|-------|--|
| 研究背景 | This gene likely encodes a member of the secernin family of proteins. A similar protein in rat functions in regulation of exocytosis in mast cells. Alternatively spliced transcript variants have been described. |
| 基因ID | 9805 |
| 基因名 | SCRN1 |
| Swiss | Q12765 (https://www.uniprot.org/uniprotkb/Q12765/entry) |
| 别名 | SES1,SCRN1 Rabbit pAb,SCRN1,KIAA0193 |

产品验证



Western blot analysis of SCRN1 expressed in Mouse brain,Rat brain,A-549 using SCRN1 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blo cking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)