

SCRN1 Rabbit pAb

货号: **AYP22867**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	46kDa
实测分子量	50KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,PC-3,U-251MG,Mouse brain,Rat brain
细胞定位	cytoplasm,nucleus,presynapse
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 195-414 of human SCR N1 (NP_055581.3).
序列	TKMDAEHPELRSYAQSQGWWTGEGEFNFSEVFPVEDHDCGAGKDSLEKQEESITVQTMNMLRDKASGVCIDSEFF LTTASGVSVLPQNRSSPCIHIFTGTPDPSRSIFKPFIFVDDVKLVPKTQSPCFGDDDDPAKKEPRFQEKPDRRHELYKAHEW ARAIIESDQEQGRKLRSTMLELEKQGLEAMEEILTSSEPLDPAEVGDLFYDCVDTEIKFFK

靶点信息

研究背景	This gene likely encodes a member of the secernin family of proteins. A similar protein in rat functions in regulation of exocytosis in mast cells. Alternatively spliced transcript variants have been described.
基因ID	9805

基因名	SCRN1
Swiss	Q12765
别名	SES1

产品验证



Western blot analysis of SCRN1 expressed in Mouse brain,Rat brain,A-549 using SCRN1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane . Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn