

— ABLYBIO, Help Your Research



CD31/PECAM1 Rabbit pAb

货号: **AYP22870**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	76kDa
实测分子量	130KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Mouse spleen,Mouse heart,Rat heart
细胞定位	cell surface,cytoplasm,external side of plasma membrane,extracellular space,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 611-678 of rat CD31/PECAM1 (NP_113779.1).
------	---

靶点信息

研究背景	cell adhesion molecule, also intracellular signaling molecule; may mediate leukocyte adhesion and migration, angiogenesis, and thrombosis
------	---

基因ID	29583
基因名	Pecam1
Swiss	Q3SWT0 (https://www.uniprot.org/uniprotkb/Q3SWT0/entry)
别名	CD31,Pecam,PECAM1,CD31/PECAM1 Rabbit pAb

产品验证



Western blot analysis of CD31/PECAM1 expressed in Mouse lung,Rat heart,Mouse spleen using CD31/PECAM1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)