

— ABLYBIO, Help Your Research



OsGSK5 Rabbit pAb

货号: **AYP23052**

产品信息

反应	Oryza sativa
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	50kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Spikes,milk grain
细胞定位	cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant protein of Oryza sativa OsGSK5.
------	---

靶点信息

研究背景	OsAAP6 greatly enhances root absorption of a range of amino acids and has effects on the distribution of various amino acids.
------	---

基因ID	4334738
基因名	
Swiss	Q6AVQ3 (https://www.uniprot.org/uniprotkb/Q6AVQ3/entry)
别名	qTGW3,qGL6,OsSK41,OsGSK5,qGL3.3,LOC4334738,OsGSK5 Rabbit pAb

产品验证



Western blot analysis of OsGSK5 expressed in Spikes,milk grain using OsGSK5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)