

EV71 3D Rabbit pAb

货号: **AYP23065**

产品信息

反应	
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	
实测分子量	51KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant protein of Enterovirus EV71 3D
序列	

靶点信息

研究背景	Enterovirus A 71 (EV71) is a neurotropic virus that may lead to acute flaccid paralysis, encephalitis, cardiopulmonary failure or even death. No vaccine and defensive drug controlling EV71 is currently available, novel and efficient antiviral drug or vaccine is therefore urgently needed. 3Dpol (RNA-dependent RNA polymerase (RdRp)) has been an important target for anti-EV71 drug development.
基因ID	
基因名	

Swiss	
别名	

产品验证



Western blot analysis of EV71 3D expressed in HeLa using EV71 3D Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn