

FAM3A Rabbit pAb

货号: **AYP23066**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	21kDa/23kDa/25kDa
实测分子量	21kDa/23kDa/25kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 34-230 of human FAM3A (NP_068578.2).
序列	FPRIQQLFTSPESSVTAAPRARKYKCGLPQPCPEEHLAFRVVSGAANVIGPKICLEDKMLMSSVKDENVGRGLNIALVNGVS GELIEARAFDMWAGDVNDLLKFIRPLHEGTLVVFASYDDPATKMNEETRKLFSSELGSRNAKELAFRDSWVFGAKGVQN KSPFEQHVKNKSHSNKYEGWPEALEMEGCIPRRSTAS

靶点信息

研究背景	This gene encodes a cytokine-like protein. The expression of this gene may be regulated by peroxisome proliferator-activated receptor gamma, and the encoded protein may be involved in the regulation of glucose and lipid metabolism. Alternative splicing results in multiple transcript variants.
基因ID	60343

基因名	FAM3A
Swiss	P98173
别名	FAM3A;2.19;DLD;DXS560S;XAP-7

产品验证



Western blot analysis of FAM3A expressed in human liver cancer using FAM3A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn