

— ABLYBIO, Help Your Research



OsSRT1 Rabbit pAb

货号: **AYP23093**

产品信息

反应	Oryza sativa
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	54kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	leaves
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant protein of Oryza sativa OsSRT1.
------	---

靶点信息

研究背景	OsSRT1 is a NAD ⁺ -dependent histone deacetylase, closely related to the human SIRT6 that plays key role in genome stability and metabolic homeostasis.
------	--

基因ID	4335342
基因名	SRT1
Swiss	Q7XWV4 (https://www.uniprot.org/uniprotkb/Q7XWV4/entry)
别名	OsSRT1,OsSRT1 Rabbit pAb,SRT1,Regulatory protein SIR2 homolog 1,SIR2-like protein 1,SRT701,LOC4335342

产品验证



Western blot analysis of OsSRT1 expressed in leaves using OsSRT1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)