

PACSIN2 Rabbit pAb

货号: AYP23500

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	51kDa/55kDa
实测分子量	51kDa/55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell membrane,Cell projection,Cytoplasm,Cytoplasmic side,Cytoplasmic vesicle membrane,Early endosome,Membrane,Peripheral membrane protein,Recycling endosome membrane,caveola,cytoskeleton,ruffle membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 240-410 of human PACSIN2 (NP_009160.2).
序列	RLRFFREVLLLEVQKHLDSLNVAGYKAIYHDLEQSIRAADAVEDLRWFRANHGPGMAMNWPQFEIEWSADLNRTLRSREKK KATDGVTLTGINQTGDQSLPSKPSSTLNVPSNPAQSAQSQSSYNPFEDDDTGSTVSEKDDTKAKNVSSYEKTQSYPTD WSDDESNNPFSST

靶点信息

研究背景	This gene is a member of the protein kinase C and casein kinase substrate in neurons family. The encoded protein is involved in linking the actin cytoskeleton with vesicle formation by regulating tubulin polymerization. Alternative splicing results in multiple transcript variants.
------	---

基因ID	11252
基因名	PACSIN2
Swiss	Q9UNF0
别名	PACSIN2; SDPII; protein kinase C and casein kinase substrate in neurons 2

产品验证



Western blot analysis of PACSIN2 expressed in human breast cancer using PACSIN2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn