

— ABLYBIO, Help Your Research



SIRT3 Rabbit pAb

货号: **AYP23527**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC
推荐浓度	IHC: 1:50 - 1:200
理论分子量	28kDa/43kDa
实测分子量	28kDa/43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HEK-293 cells, HeLa cells, HepG2 cells
细胞定位	Mitochondrion matrix
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 300-399 of human SIRT3 (NP_03637 1.1).
------	---

靶点信息

研究背景	This gene encodes a member of the sirtuin family of proteins, homologs to the yeast Sir2 protein. Members of the sirtuin family are characterized by a sirtuin core domain and grouped into four classes. The functions of human sirtuins have not yet been determined; however, yeast sirtuin proteins are known to regulate epigenetic gene silencing and suppress recombination of rDNA. Studies suggest that the human sirtuins may function as intracellular regulatory proteins with mono-ADP-ribosyltransferase activity. The protein encoded by this gene is included in class I of the sirtuin family. Two alternatively spliced transcript variants that encode different proteins have been described for this gene.
基因ID	23410
基因名	SIRT3
Swiss	Q9NTG7 (https://www.uniprot.org/uniprotkb/Q9NTG7/entry)
别名	SIRT3,SIR2L3,sirtuin 3,SIRT3 Rabbit pAb,NAD-dependent protein deacetylase sirtuin-3,Regulatory protein SIR2 homolog 3,SIR2-like protein 3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)