

ANKRD52 Rabbit pAb

货号: **AYP23539**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	115kDa
实测分子量	115kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 450-660 of human ANKRD52 (NP_775866.2).
序列	RRRDKFGRTPLHYAAANGSYQCAVTLVTAGAGVNEADCKGCSPLHYAAASD TYRRAEPHTPSSHDAEED EPLKESRRKEA FFCLEFLLDNGADPSLRDRQGYTAVHYAAAYGNRQNLLELLEMSFNCL EDVESTIPVSPLHLAAYNGHCEALKTLAETLVNL DVRDHKGR TALFLATERGSTECEVLT AHGASALIKERKRKWTPLHAAA

靶点信息

研究背景	Putative regulatory subunit of protein phosphatase 6 (PP6) that may be involved in the recognition of phosphoprotein substrates.
基因ID	283373

基因名	ANKRD52
Swiss	Q8NB46
别名	ANKRD52; ANKRD33; ankyrin repeat domain 52

产品验证



Western blot analysis of ANKRD52 expressed in human lung cancer using ANKRD52 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn