

# NRF2 Rabbit pAb

货号: **AYP23578**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	65kDa/66kDa/67kDa
实测分子量	97-100KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C6
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant protein of human NRF2.
序列	

## 靶点信息

研究背景	This gene encodes a transcription factor which is a member of a small family of basic leucine zipper (bZIP ) proteins. The encoded transcription factor regulates genes which contain antioxidant response elements (ARE) in their promoters; many of these genes encode proteins involved in response to injury and inflammation which includes the production of free radicals. Multiple transcript variants encoding different isoforms have been characterized for this gene.
基因ID	4780

基因名	NFE2L2
Swiss	Q16236
别名	NFE2L2; HEBP1; NRF2; nuclear factor; erythroid 2 like 2

## 产品验证



Western blot analysis of NRF2 expressed in C6 using NRF2 Rabbit pAb at 1:1000. Secondary antibody: HR P Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)