

HNRNPA2B1 Rabbit pAb

货号: **AYP23717**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	36kDa/37kDa
实测分子量	36kDa/37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A431 cells, Jurkat cells, Raji cells, A549 cells, HL-60 cells, HeLa cells, K-562 cells, HSC-T6 cells
细胞定位	Cytoplasm,Nucleus,nucleoplasm
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human HNRNP A2B1 (NP_112533.1).
序列	MEKTLETVPLERKKREKEQFRKLFIGGLSFETTEESLRNYEQWGKLTDCVVMRDPASKRSRGFGFVTFSSMAEVDAAAMA ARPHSIDGRVVEPKRAVAREESGKPGAHVTVKKLFVGGIKEDTEEHHLRDYFEEYGKIDTIEITDRQSGKKRGGFGFVTFDD HDPVDKIVLQKYHTINGH

靶点信息

研究背景	This gene belongs to the A/B subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has two repeats of quasi-RRM domains that bind to RNAs. This gene has been described to generate two alternatively spliced transcript variants which encode different isoforms.
基因ID	3181
基因名	HNRNPA2B1
Swiss	P22626
别名	HNRNPA2B1; HNRNPA2; HNRNPB1; HNRPA2; HNRPA2B1; HNRPB1; IBMPFD2; RNPA2; SNRPB1; heterogeneous nuclear ribonucleoprotein A2/B1

产品验证



Western blot analysis of HNRNPA2B1 expressed in A431 cells, Jurkat cells, Raji cells using HNRNPA2B1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn