

# NLRP3 Rabbit pAb

货号: **AYP23758**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:1000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	
实测分子量	110KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	RAW264.7
细胞定位	cytoplasm,cytosol,endoplasmic reticulum,extracellular region,inflammasome complex,NLRP3 inflammasome complex,nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-360 of mouse NLRP3 (NP_665826.1).
序列	MTSVRCKLAQYLEDLEDVDLKKFKMHLEDYPPEKGCIPVPRGQMEKADHLDLATLMIDFNGEEKAWAMAVWIFAAINRRDLWEKAKKDQPEWNDTCTSHSSMVCQEDSLEEEWMLLGYLSRISICKKKKDYCKMYRRHVRSRFYSIKDRNARLGESVDLNSRYTQLQLVKEHPSKQEREHELLTIGRTKMRDSPMSSLKLELLFEPEDGHSEPVHTVVFQGAAGIGKTLARKIMLDWALGKLFKDKFDYLFHICREVSLRTPRSLADLIVSCWPDNPPVCKILRKPSRILFLMDGFDELQGAFFEDEHIGEVCTDWQKAVRGDILLSSLRKKLLPKASLLITTRPVALEKLQHLLDH

## 靶点信息

研究背景	This gene encodes a pyrin-like protein containing a pyrin domain, a nucleotide-binding site (NBS) domain, and a leucine-rich repeat (LRR) motif. This protein interacts with the apoptosis-associated speck-like protein PYCARD/ASC, which contains a caspase recruitment domain, and is a member of the NLRP3 inflammasome complex. This complex functions as an upstream activator of NF-kappaB signaling, and it plays a role in the regulation of inflammation, the immune response, and apoptosis. The SARS-CoV 3a protein, a transmembrane pore-forming viroporin, has been shown to activate the NLRP3 inflammasome via the formation of ion channels in macrophages. Mutations in this gene are associated with familial cold autoinflammatory syndrome (FCAS), Muckle-Wells syndrome (MWS), chronic infantile neurological cutaneous and articular (CINCA) syndrome, neonatal-onset multisystem inflammatory disease (NOMID), keratoendotheliitis fugax hereditaria, and deafness, autosomal dominant 34, with or without inflammation. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. Alternative 5' UTR structures are suggested by available data; however, insufficient evidence is available to determine if all of the represented 5' UTR splice patterns are biologically valid.
基因ID	216799
基因名	
Swiss	Q8R4B8
别名	ALI; AVP; FCU; MWS; FCAS; KEFH; CIAS1; FCAS1; NALP3; C1orf7; CLR1.1; DFNA34; PYPAF1; AGTAVPRL

## 产品验证



Western blot analysis of NLRP3 expressed in RAW264.7 using NLRP3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)