

NPHS1 Rabbit pAb

货号: **AYP23821**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	130kDa/134kDa
实测分子量	200KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549
细胞定位	Cell membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1019-1241 of human NP HS1 (NP_004637.1).
序列	LGDSGLADKGTQLPITTPGLHQPSGEPEDQLPTEPPSGPSGLPLLPLVLFALGGLLLSNASCVGGVWLQRRRLRLAEGISE KTEAGSEEDRVRNEYEESQWTGERDTSSTVSTTEAEPYRSLRDFSPQLPPTQEEVSYSRGFTGEDEDMAPFGHLYDE VERTYPPSGAWGPLYDEVQMGPWDLHWPEDTYQDPRGIYDQVAGDLDTLEPDSLPELRGHLV

靶点信息

研究背景	This gene encodes a member of the immunoglobulin family of cell adhesion molecules that functions in the glomerular filtration barrier in the kidney. The gene is primarily expressed in renal tissues, and the protein is a type-1 transmembrane protein found at the slit diaphragm of glomerular podocytes. The slit diaphragm is thought to function as an ultrafilter to exclude albumin and other plasma macromolecules in the formation of urine. Mutations in this gene result in Finnish-type congenital nephrosis 1, characterized by severe proteinuria and loss of the slit diaphragm and foot processes.
------	---

基因ID	4868
基因名	NPHS1
Swiss	O60500
别名	NPHS1;CNF;NPHN;nephrin

产品验证



Western blot analysis of NPHS1 expressed in A-549 using NPHS1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn