

— ABLYBIO, Help Your Research



# Aim2 Rabbit pAb

货号: **AYP24524**

## 产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	40kDa
实测分子量	39KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat spleen
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant protein of mouse Aim2.
------	------------------------------------

## 靶点信息

研究背景	Enables double-stranded DNA binding activity and identical protein binding activity. Involved in positive regulation of macromolecule metabolic process; pyroptosis; and regulation of defense response. Acts upstream of or within cellular response to interferon-beta. Located in cytoplasm. Part of AIM2 inflammasome complex. Orthologous to human AIM2 (absent in melanoma 2).
基因ID	383619
基因名	Aim2
Swiss	Q91VJ1
别名	Gm1313; Ifi210

## 产品验证



Western blot analysis of Aim2 expressed in Rat spleen using Aim2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)