

— ABLYBIO, Help Your Research



CGPR Rabbit pAb

货号: **AYP24795**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	13kDa
实测分子量	16KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse thymus
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant protein of Human CGPR.
------	------------------------------------

靶点信息

研究背景	Predicted to enable calcitonin receptor binding activity. Predicted to be involved in adenylate cyclase-activating G protein-coupled receptor signaling pathway and regulation of cytosolic calcium ion concentration. Predicted to be located in extracellular region. Predicted to be active in extracellular space.
------	--

基因ID	797
基因名	CALCB
Swiss	P10092 (https://www.uniprot.org/uniprotkb/P10092/entry)
别名	CALC2, CGRP2, CGRP-II, CGPR Rabbit pAb, CALCB, Beta-type CGRP, Calcitonin gene-related peptide II

产品验证



Western blot analysis of CGPR expressed in Mouse thymus using CGPR Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)