

CGPR Rabbit pAb

货号: **AYP24795**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	13kDa
实测分子量	16KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse thymus
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant protein of Human CGPR.
序列	

靶点信息

研究背景	Predicted to enable calcitonin receptor binding activity. Predicted to be involved in adenylate cyclase-activating G protein-coupled receptor signaling pathway and regulation of cytosolic calcium ion concentration . Predicted to be located in extracellular region. Predicted to be active in extracellular space.
基因ID	797
基因名	CALCB

Swiss	P10092
别名	CALC2; CGRP2; CGRP-II

产品验证



Western blot analysis of CGPR expressed in Mouse thymus using CGPR Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn