

# BRD7 Rabbit pAb

货号: **AYP24823**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	74kDa
实测分子量	75KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,U266,Mouse brain,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant protein of human BRD7.
序列	MGKHKHKKHSDKHL YEEYVEKPLKLVLVGGNEVTELTSGSSGHDSSLFEDKNDHKHKDRKRKRKKGKQIPGEEKG RKR RRVKEDKKKRDRD

## 靶点信息

研究背景	This gene encodes a protein which is a member of the bromodomain-containing protein family. The product of this gene has been identified as a component of one form of the SWI/SNF chromatin remodeling complex, and as a protein which interacts with p53 and is required for p53-dependent oncogene-induced senescence which prevents tumor growth. Pseudogenes have been described on chromosomes 2, 3, 6, 13 and 14. Alternative splicing results in multiple transcript variants.
基因ID	29117

基因名	BRD7
Swiss	Q9NPI1
别名	BRD7;BP75;CELTIX1;NAG4

## 产品验证



Western blot analysis of BRD7 expressed in HepG2,U266,Mouse brain using BRD7 Rabbit pAb at 1:1000 . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)