

# FOXO3A Rabbit pAb

货号: **AYP25138**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus , Homo sapiens , Rattus norvegicus , Danio rerio <b>Co-IP:</b> Rattus norvegicus <b>IP:</b> Rattus norvegicus
应用	<a href="#">WB</a> <a href="#">IHC</a> <a href="#">IF/ICC</a>
推荐浓度	<b>WB:</b> 1:500 - 1:1000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	71kDa
实测分子量	80-90KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,293T,PC-3
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of mouse FOXO3 A (NP_062714.1).
序列	MAEAPASPVPPLSPLEVELDPEFEPQSRPRRSTWPLQRPELQASPAKPSGETAADSMIPEEDDDDEDGGRASSAMVIGGGVSSTLGSGLLLED SAMLLAPGGQDLGSGPASAAGALSGGTPTQLQPQPPLQPQPGAAGGSGQPRKCSSRRNAW

## 靶点信息

研究背景	This gene belongs to the forkhead family of transcription factors which are characterized by a distinct fork head domain. This gene likely functions as a trigger for apoptosis through expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding the same protein have been observed.
基因ID	56484
基因名	
Swiss	Q9WVH4
别名	FOXO3;AF6q21;FKHRL1;FKHRL1P2;FOXO2;FOXO3A

## 产品验证



Western blot analysis of FOXO3A expressed in HeLa, Jurkat, 293T using FOXO3A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)